

# CARDIOVIT AT-2

Exceed your Expectations,  
not your Budget!



## CARDIOVIT AT-2 – Why compromise Performance and Reliability when you can have it all in one affordable ECG-System?


The AT-2 is a low cost ECG system with performance features previously only available in high priced systems. Using the latest design and manufacturing technology, the AT-2 provides A4 reports within seconds at the touch of a button. Optional is an automatic measurement and interpretation program as well as a trolley.

- Full-size A4 reports to fit patient charts
- Powerful battery, which lasts for 400 recordings or 6 hours continuous use
- Built-in thermal printer
- Up to 20 different printout formats
- Realtime printout of 6 channels with 5, 10, 25 oder 50 mm/s
- Various recording rates
- Simultaneous 12 lead ECG acquisition
- Printouts in less than 20 seconds for immediate analysis
- Auto printout of all **12 channels** on 1 or 2 pages
- Copy function for unlimited number of copies
- AC/DC power supply with built-in rechargeable battery, built-in charger unit
- Portable – weighs only 4.8 kg
- CE labelled
- Meets or exceeds U.S. and international quality standards

### Options:

- *Measurement and interpretation software*
- *Cart*

Seeing is believing... call now for a no-obligation demonstration!



# SCHILLER

www.schiller.ch

The Art of Diagnostics