


Select the instruments you need to perform a complete diagnostic examination:



767 Wall Transformer

Handle-based rheostats and automatic on/off sensors make operation simple.



PanOptic™ Ophthalmoscope

5x larger field of view of the fundus vs. standard ophthalmoscopes in an undilated eye (25° vs. 5°).



3.5 V Coaxial Ophthalmoscope

Multiple aperture/filter options provide long-lasting, reliable performance for general and specialist examinations.



MacroView™ Otoscope

A nearly complete view of the tympanic membrane with approximately twice the field of view and 30% greater magnification than a traditional otoscope.



Pneumatic Otoscope

Halogen HPX™ lamp that provides 30% more light output for true tissue color and consistent, long-lasting illumination.



767 Wall Aneroid Sphygmomanometer

Easy-to-read, jeweled-movement, mercury-free manometer is certified to be accurate within one percent of scale (± 3 mmHg).



KleenSpec® Disposable Specula Dispenser

Offers optimum storage and organization of disposable otoscope tips and other supplies. Comes filled with specula.



SureTemp® Plus Electronic Thermometer

Provides accurate oral temperature readings in about four to six seconds, 10- to 13-second pediatric axillary and rectal temperature readings and 12- to 15-second adult axillary temperature readings.



Spot Vital Signs®/Spot Vital Signs® LXi Multiparameter Spot Check Devices

Fast, accurate, and efficient way to take vitals. Allows for electronic transfer of data to EHR.

© 2011 Welch Allyn MC8107

The Welch Allyn 767 Integrated Diagnostic System

Recommended Systems:



76791-MX

Includes 767 Transformer (#76710), Coaxial Ophthalmoscope (#11720), MacroView Otoscope (#23810), 767 Aneroid (#7670-01), KleenSpec Dispenser—small (#52401), and SureTemp Plus Thermometer (#01690-300)



76791-2MPX

Includes 767 Transformer (#76710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), 767 Aneroid (#7670-01), KleenSpec Dispenser—small (#52401), and SureTemp Plus Thermometer (#01690-300)



76792-M

Includes 767 Transformer (#76710), Coaxial Ophthalmoscope (#11720), MacroView Otoscope (#23810), 767 Aneroid (#7670-01), and KleenSpec Dispenser—small (#52401)



76794-2MPLXI

Includes 767 Transformer (#76710), PanOptic Ophthalmoscope (#11820), MacroView Otoscope (#23810), KleenSpec Dispenser—small (#52401), and Wall Mount for Spot Vital Signs LXi (#4500-62). Spot Vital Signs/Spot Vital Signs LXi ordered separately.

WelchAllyn®

Advancing Frontline Care™

All the instruments needed to perform a complete diagnostic examination.

	WALL BOARD	767 WALL TRANSFORMER	PANOPTIC™ OPHTHALMOSCOPES		COAXIAL OPHTHALMOSCOPE	DIAGNOSTIC MACROVIEW™ OTOSCOPIES		PNEUMATIC OTOSCOPE	767 ANEROID	KLENSPEC® DISPENSERS		SURETEMP® PLUS THERMOMETER
Part #		76710	11810	11820	11720	23810	23820	20200	7670-00	52401 (small)	52101 (large)	01690-300
76791-MX	•	•			•	•			•	•		•
76791-1MPX	•	•	•			•			•	•		•
76791-2MPX	•	•		•		•			•	•		•
76792-M	•	•			•	•			•	•		
76792-MNOBP	•	•			•	•				•		
76792-1MP	•	•	•			•			•	•		
76792-2MP	•	•		•		•			•	•		
76793-MX	•	•			•	•			•	•		•
76794-M	•	•			•	•			•	•		
76795-MX	•	•			•		•		•	•		•
76795-1MPX	•	•	•				•		•	•		•
76795-2MPX	•	•		•			•		•	•		•
76796-M	•	•			•		•		•	•		
76796-1MP	•	•	•				•		•	•		
76796-2MP	•	•		•			•		•	•		
76797-X	•	•			•			•	•		•	•
76797-2PX	•	•		•				•	•		•	•
79798-2P	•	•		•				•	•		•	
76794-MSPOT*	•	•			•	•				•		
76794-MLXI*	•	•			•	•				•		
76794-2MPSPOT*	•	•		•		•				•		
76794-2MPLXI*	•	•		•		•				•		

*Spot Vital Signs®/Spot Vital Signs® LXi automated vital signs devices must be ordered separately. Most popular systems listed. See your Welch Allyn representative for a complete listing.