

# EASY, TOUCH VITALS BY MIDMARK



IQ  
vitals

MIDMARK®

Midmark IQvitals™ is an affordable, Windows® based, touchscreen device that provides automated blood pressure, temperature, SpO<sub>2</sub> and pulse rate readings with a touch of the screen. An alphanumeric keyboard allows for fast and easy entry of patient data, and the intuitive design increases productivity and simplifies staff training. IQvitals™ can be used as a stand-alone unit or will connect easily to a computer using IQmanager™ software.

# IQvitals™ Technology by Midmark

Use as a standalone unit or connect to a computer using IQmanager™ software

## TECHNOLOGICALLY ADVANCED, USER DESIGNED VITAL SIGNS

### Efficiency

- Automated, all-in-one vital signs device saves time.
- Seamless connection to IQmanager™ software reduces transcription errors.
- Lightweight, portable device weighs just 3.9 lbs.
- Long battery life allows vitals to be captured continuously and quickly.
- Connects to select digital scales – use either digital or manual option to enter weight.
- Easily enter and view all vital signs data in one location, saving time.
- Manual vitals entry options for respiration rate, height, weight and pain score; calculates BMI automatically.

### Technology

- Windows® based platform provides flexible, easy-to-upgrade software for the life of the product.
- Simplified touchscreen display for effortless, one-touch operation.
- Alphanumeric keyboard allows for quick and easy patient data entry.
- Intuitive user interface simplifies staff training.
- Meets or exceeds all standards for blood pressure, temperature, pulse and SpO<sub>2</sub>.
- Built with hospital grade components to stand up to the demands of a busy practice.

### Adaptability

- USB and serial connectivity options are standard.
- 100-reading memory for easy review, printing or download.
- Accessories available for patient testing, from pediatric to adult.
- Mobile cart, table top or wall mounting options available.
- External printer available.



Flexible mounting options to meet your needs.

Automated, all-in-one vital signs device saves time and reduces transcription errors.



Full-color, touchscreen display is quick and easy to operate.



Lightweight and portable at just 3.9 pounds with table-top, wall or mobile cart mounting options.



Connected flexibility with USB and serial ports standard.

## Make your practice more efficient with an IQ Digital Diagnostic System

Workflow efficiency and improvements in patient care are just a few of the benefits that can result from implementing an IQ Digital Diagnostic System. Call us today and we'll put together a customized suite of solutions to help grow as your practice grows.

Diagnostic Products									
	ECG	Holter	Stress	Spirometer	Vital Signs	EHR/EMR Connectivity	Single Platform (IQmanager™)	Stand-alone	Thin - Client (IQpath™)
IQecg™	■					■	■		■
IQflex™	■							■	
IQclassic™	■							■	
IQholter™		■				■	■		
IQholter™ EX		■				■	■		
IQholter™ EP		■				■	■		
IQstress™			■			■	■		
IQspiro™				■		■	■		■
IQvitals™					■		■	■	
IQvitals™ PC					■		■		

### System Requirements

- Windows® Vista Home, Windows® XP Pro, Windows® 7 Home Premium, 32-bit and 64-bit
- Intel® Core™ 2 Duo Processor E4300 or 64-bit processor or faster
- 2 GB of RAM
- 2 GB available hard disk space
- USB or Serial Port
- Mouse, keyboard, and high resolution display monitor

Please see the IQvitals™ Operation Manual or visit [midmark.com](http://midmark.com) for more detailed system requirements

### Ordering Information

- 4-000-0500 IQvitals™ with Blood Pressure and Alaris® Turbo Temp® Thermometer
- 4-000-0510 IQvitals™ with Blood Pressure, Alaris® Turbo Temp® Thermometer and SpO<sub>2</sub>
- 1-100-1605 IQvitals™ Thermal Printer
- 1-100-1603 Fairbanks® TeleWeigh™ Digital Scale
- 3-004-2000 IQvitals™ Mobile Cart

*Alaris and Turbo Temp are registered trademarks of Cardinal Health, Inc.  
Windows is a registered trademark of Microsoft Corporation.  
Intel Core is a trademark of Intel Corporation in the U.S. and/or other countries.  
Fairbanks and TeleWeigh are trademarks of Fairbanks Scales, Inc.*

For information, pricing or to place an order, contact us at:  
1-800-624-8950 or  
1-310-516-6050 or  
[midmark.com](http://midmark.com)

Midmark Diagnostics Group is an ISO 13485 Certified Company.

Midmark Diagnostics Group  
1125 W. 190th Street, Gardena, CA 90248  
[midmark.com](http://midmark.com)  
1-800-624-8950

© 2010 Midmark Corporation  
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286  
Products subject to improvement changes without notice  
Litho in U.S.A. 007-0836-00 Rev. A1 (8/10)



## RETHINK THE CLINICAL SPACE

At the core of Midmark is a new way to think about caregiver/patient interaction and the delivery of efficient patient care. We are dedicated to working together with you to redefine the future of the clinical space.

[midmarkclinicalsolutions.com](http://midmarkclinicalsolutions.com)

